



# AFIAS BNP

## INTENDED USE

**AFIAS BNP** is a fluorescence Immunoassay (FIA) for the quantitative determination of cardiac B-type Natriuretic Peptide(BNP) level in human whole blood/plasma. It is useful as an aid in management and monitoring of Heart Failure(HF).  
 For *in vitro* diagnostic use only.

## INTRODUCTION

B-type natriuretic peptide (BNP) is a cardiac neurohormone that derives from the precursor pre-proBNP, containing 134 amino acids and including a signal peptide of 26 amino acids. proBNP, produced by cleavage of the signal peptide, is further split into BNP, which is considered to be the biologically active hormone, and an inactive amino terminal fragment, NT-proBNP. Blood measurements of BNP have been shown to be useful in differentiating dyspnoea caused by congestive heart failure (CHF) from dyspnoea related to other causes.

**AFIAS BNP** Assay quantitatively measures the concentration of BNP, providing a physician with an accurate tool to assess prognosis in patients with chronic heart failure.

## PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes, and migrate onto nitrocellulose matrix to be captured by the other immobilized-antibodies on the test strip.

More antigens in sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by the instrument for AFIAS tests to show BNP concentration in the sample.

## COMPONENTS

- AFIAS BNP consists of 'cartridges'.
- Each sealed aluminum pouch contains two cartridges.
- Each cartridge packaged in an aluminum pouch has three components including a cartridge part, detector part and diluent part.
- The cartridge part contains the membrane called a test strip which has streptavidin at the test line and chicken IgY at the control line.
- The detector part has 2 granules containing anti-BNP fluorescence conjugate, anti-BNP biotin conjugate, anti-chicken IgY-fluorescence conjugate, bovine serum albumin, bromophenol blue and sodium azide as a preservative in Tris buffer.
- The diluent part contains tween 20, sodium chloride in sodium phosphate buffer.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow the instructions and procedures described in this 'Instructions for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges, c-tips, or tips. A cartridge should be used for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use a cartridge, if the pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- If test components and/or sample are stored in refrigerator, then allow them to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges, c-tips and pipette tips should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- The cartridge contains sodium azide (NaN<sub>3</sub>), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.
- No Biotin interference was observed in **AFIAS BNP** when biotin concentration in the sample was below 5 ng/ml. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
- AFIAS BNP** will provide accurate and reliable results subject to the below conditions.

- **AFIAS BNP** should be used only in conjunction with the instrument for AFIAS tests.

Recommended anticoagulant
K <sub>3</sub> EDTA

### [30 µL C-tip]

- Wear disposable gloves and protective equipment for safety.
- Do not reuse the c-tips as they are disposable.
- Check the surface for damage or contamination.
- Careful when collecting sample to prevent air bubbles from forming in the c-tip.
- Careful not to get blood on the surface of the c-tip. Wipe the surface of the c-tip with tissue.
- Avoid placing them under direct sunlight and keep them in dry places.
- Unused c-tips must be stored in a sealed zipper bag
- Sample collection tool and sample container are infectious and therefore must be handled carefully and discarded by an appropriate method in accordance with relevant local

regulations. .

## LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes the antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

## STORAGE AND STABILITY

Component	Storage condition		
	Storage Temperature	Shelf life	Note
Cartridge	2 – 30 °C	20 months	Unopened
	2 – 30 °C	1 month	Resealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

## MATERIALS SUPPLIED

REF SMFP-83 Components of AFIAS BNP	
Cartridge Box Contains	
- Cartridge	24
- Pipette tip (zipper bag)	24
- Spare cartridge zipper bag	1
- ID chip	1
- Instruction for use	1

## MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS BNP**.

Please contact our sales division for more information.

AFIAS-1	REF FPRR019
AFIAS-6	REF FPRR020
AFIAS-10	REF FPRR038
AFIAS-3	REF FPRR040
C-tip (Zipper bag, 30 µL)	REF CFPO-199
Boditech BNP Control	REF CFPO-304
Boditech BNP Calibrator	REF CFPO-305

## SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS BNP** is human whole blood/plasma.

- It is recommended to test the plasma sample within 8 hours at room temperature and within 24 hours at 2-8 °C after collection.
- The plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood. If longer storage is required, e.g. if the test could not be performed within 8 hours at room temperature and within 24 hours at 2-8 °C, plasma should be immediately frozen below -20 °C. The freezing storage of sample up to 12 months does not affect the quality of results.
- However, the whole blood sample should not be kept in a freezer in any case.
- After thawing, gently shake the sample to mix.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.
- Collection of whole blood sample using C-tip
  - Hold the C-tip horizontally and touch the surface of the blood with the tip of the C-tip.
  - Capillary action will automatically draw the blood sample to C-tip and stop.
  - Wipe off any excess blood around the tip.
  - Double-check if whole blood is filled accurately in the C-tip and AFIAS reader is ready for a test on the 'C-tip mode'.

## TEST SETUP

- Check the components of the **AFIAS BNP** as described below. : Cartridges, pipette tips, C-tips, an ID chip, a spare cartridge zipper bag and an instructions for use.
- Ensure that the lot number of the cartridge matches that of the ID chip.
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the 'ID chip port'.
- ※ **Please refer to the instrument for AFIAS tests operation manual for complete information and operating instructions.**

## TEST PROCEDURE

### ▶ AFIAS-6, AFIAS-1, AFIAS-3

#### General mode

- Insert a cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Select the 'General mode' in the instrument for AFIAS tests.
- Take 100 µL of the sample (whole blood /plasma/control) using a pipette and dispense it into the sample well of the cartridge.
- Tap the 'Start' button on the screen.
- The test result will be displayed on the screen after 12 minutes.

#### C-tip mode

- Insert a cartridge into the cartridge holder.

- 2) Take 30 µL of whole blood using a C-tip.
- 3) Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- 4) Select the 'C-tip mode' in the instrument for AFIAS tests.
- 5) Tap the 'Start' button on the screen.
- 6) The test result will be displayed on the screen after 12 minutes.

▶ **AFIAS-10**

**Normal mode**

- 1) Insert a cartridge into the cartridge holder.
- 2) Insert a tip into the tip hole of the cartridge.
- 3) Tap the 'Load' button of the bay that holds the cartridge with the tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 4) Insert the sample tube into the tube rack.
- 5) Insert the tube rack into the loading part of the sampling station.
- 6) Tap the 'Start' button on the screen.
- 7) The test result will be displayed on the screen after 12 minutes.

**Emergency mode – General tip**

- 1) The test procedure is same with the 'Normal mode 1) – 3)'.  
2) Convert the 'Emergency mode' in AFIAS-10.
- 3) Select the tip type (general tip) on the screen.
- 4) Select the sample type (whole blood /plasma) on the screen.
- 5) Take 100 µL of the sample using a pipette and dispense it into the sample well of the cartridge.
- 6) Tap the 'Start' button on the screen.
- 7) The test result will be displayed on the screen after 12 minutes.

**Emergency mode – C-tip**

- 1) Insert a cartridge into the cartridge holder.
- 2) Take 30 µL of whole blood using a C-tip.
- 3) Insert the C-tip with sample into the tip hole of the cartridge.
- 4) Tap the 'Load' button of the bay that holds the cartridge with a tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 5) Convert the 'Emergency mode' in AFIAS-10.
- 6) Select the tip type (C-tip) on the screen.
- 7) Tap the 'Start' button on the screen.
- 8) The test result will be displayed on the screen after 12 minutes.

※ Note: Refer to the instrument for AFIAS test Operation Manual to select a sample type.

**INTERPRETATION OF TEST RESULT**

- The instrument for AFIAS tests calculate the test result automatically and displays BNP concentration of the test sample in terms of pg/mL.
- Reference value: 100 pg/mL.
- Working range: 5 – 5,000 pg/mL.

**QUALITY CONTROL**

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.

- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **AFIAS BNP**. For more information regarding obtaining the control materials, contact [Boditech Med Inc.'s Sales Division for assistance](#). (Please refer to the instruction for use of control material.)

**PERFORMANCE CHARACTERISTICS**

- **Analytical sensitivity**
  - Limit of Blank (LoB) 2.45 pg/mL
  - Limit of Detection (LoD) 3.65 pg/mL
  - Limit of Quantification (LoQ) 5.00 pg/mL

- **Analytical Specificity**
  - **Cross-reactivity**  
Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS BNP** test results did not show any significant cross-reactivity with these biomolecules.

Cross reactants	Conc.
CK-MB	1,000 ng/ml
Troponin complex	1,000 ng/ml
NT-proBNP	1,000 ng/ml
Myoglobin	200 ng/ml

- **Interference**  
Interferents listed in the following table were added to the test sample at the concentration mentioned below. **AFIAS BNP** test results did not show any significant interference with these materials.

interference materials	Conc.
L-Ascorbic acid	350 µmol/L
Bilirubin	350 µmol/L
Cholesterol	13 mmol/L
D-Glucose	1,000 mg/dL
Hemoglobin	2 g/L
Triglyceride mixture	500 mg/dL
EDTA	3.4 µmol/L

- **Precision**
  - Between Lot  
One person tested three different lots of **AFIAS BNP**, ten times at each concentration of the control standard.
  - Between person  
Three different persons tested one lot of **AFIAS BNP**, ten times at each concentration of the control standard.
  - Between day  
One person tested one lot **AFIAS BNP**, during five days, ten times at each concentration of the control standard.
  - Between site  
One person tested **AFIAS BNP** at three different sites; ten times at each concentration of the control standard.

Conc. [pg/ml]	Between lot		Between person	
	AVG	CV (%)	AVG	CV (%)
10	10.01	5.2	9.91	5.4
156	156.60	6.2	156.16	5.9
1250	1252.62	4.7	1254.46	5.7


Conc. [pg/ml]	Between day		Between site	
	AVG	CV (%)	AVG	CV (%)
10	10.04	5.1	9.97	6.0
156	156.09	5.8	158.69	5.6
1250	1235.01	6.0	1238.42	5.6

- **Accuracy**  
The accuracy was confirmed by testing 3 different lots, ten times at each concentration of the control standard.

Conc. [pg/mL]	Lot 1	Lot 2	Lot 3
10	9.99	10.32	10.08
25	23.89	25.08	25.96
50	49.10	51.62	50.19
100	102.50	102.25	99.67
250	250.96	249.39	252.12
500	501.64	476.36	489.10
1000	1014.81	993.12	1009.10
2000	2011.20	2016.87	2001.55
4500	4528.32	4465.50	4397.93

Conc. [pg/mL]	Mean	CV (%)	Recovery (%)
10	10.13	5.6	101.3
25	24.98	6.4	99.9
50	50.31	5.9	100.6
100	101.47	5.6	101.5
250	250.82	5.4	100.3
500	489.04	5.5	97.8
1000	1005.68	5.3	100.6
2000	2009.87	4.9	100.5
4500	4463.92	5.3	99.2

- **Comparability**  
BNP concentration of 100 clinical samples were quantified independently with **AFIAS BNP (AFIAS-6)** and comparator A as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and correlation coefficient (R). The regression equation and correlation coefficient are as follows.



**REFERENCES**

1. Maisel AS, Krishnaswamy PA, Nowak RM, McCord J, Hollander JE, Duc P, et al. Rapid measurement of B-type natriuretic peptide in the emergency diagnosis of heart failure. *N Engl J Med* 2002;347:161-7.
2. Mueller, T. Diagnostic accuracy of B type natriuretic peptide and amino terminal proBNP in the emergency diagnosis of heart failure. *Heart*, 2005, 91(5), 606–612.
3. Grantham, J. A., Borgeson, D. D., & Burnett, J. C. BNP: pathophysiological and potential therapeutic roles in acute congestive heart failure. *American Journal of Physiology-*

Regulatory, Integrative and Comparative Physiology, 1997, 272(4), R1077–R1083.

4. Dao, Q., Krishnaswamy, P., Kazanegra, R., Harrison, A., Amirnovin, R., Lenert, L., ... Maisel, A. S. Utility of B-type natriuretic peptide in the diagnosis of congestive heart failure in an urgent-care setting. *Journal of the American College of Cardiology*, 2001, 37(2), 379–385

**Note:** Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance; please contact:  
**Boditech Med Inc.'s Technical Services**  
Tel: +(82) 33 243-1400  
E-mail: sales@boditech.co.kr

**Boditech Med Incorporated**  
43, Geodudanji 1-gil, Dongnae-myeon,  
Chuncheon-si, Gang-won-do, 24398  
Republic of Korea  
Tel: +(82) -33-243-1400  
Fax: +(82) -33-243-9373  
www.boditech.co.kr

**Obelis s.a**  
Bd. Général Wahis 53,  
1030 Brussels, Belgium  
Tel: +(32) -2-732-59-54  
Fax: +(32) -2-732-60-03  
E-Mail: mail@obelis.net

