


Check Test Cartridge Integrity/ Kontrola neporušenosti testovací  
patrony / Kontroller, at analysebægeret ikke er beskadiget /  
Integrität der Testkartusche überprüfen/ Comprobación de  
la integridad del cartucho de análisis / Tarkista testikasetin  
käyttökelpoisuus / Vérifier l'intégrité de la cartouche de test /  
Controllo dell'integrità della cartuccia / Controleer de integriteit van  
de testcartridge / Kontrollera testkassetten integritet / Verificar  
a Integridade do Cartucho de Teste / Kontroller testkassetten  
integritet / Kontrola integrity testovacej kyvety

**EKF** | Diagnostics  
for life

EKF-diagnostic GmbH  
Ebendorfer Chausse 3  
39179 Barleben  
Germany  
T +49 39203 511 0  
F +49 39203 511 171  
info@ekfdiagnostics.com  
www.ekfdiagnostics.com

## Quo-Test A1C Test Kit

**REF** 0031


(EN)	1...10	Instructions for Use
(CZ)	11...20	Návod k použití
(DA)	21...30	Brugsanvisning
(DE)	31...40	Gebrauchsanleitung
(ES)	41...50	Instrucciones de uso
(FI)	51...60	Käyttöohjeet
(FR)	61...70	Mode d'emploi
(IT)	71...80	Istruzioni per l'uso
(NL)	81...90	Gebruiksaanwijzing
(NO)	91...100	Bruksanvisning
(PT)	101...110	Instruções de Utilização
(SV)	111...120	Bruksanvisning
(SK)	121...130	Návod na použitie

**EKF** | Diagnostics  
for life

SYMBOL	MEANING
	Consult Instructions For Use / User Manual
	In vitro Diagnostic Medical Device
	Batch or Lot Number
	Calibration Barcode
	Product Reference
	Control
	Storage Temperature
	Contents Sufficient for "Σ" Number of Tests
	Expiry Date
	Manufacturer
	Caution
	Conformity to the European Directive 98/79/EC on in vitro Diagnostic Medical Devices

References / Reference / Litteratur / Literaturhinweise / Bibliografía / Viittet / Références / Bibliografia / Referenties / Referanser / Referências / Referenser / Literatúra

1. Frantzen F, Studies on glycated haemoglobin and new analytical methods for its determination. Dr. Philos. Thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5
2. The Diabetes Control and Complications trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. N Engl J Med 1993; 329:977-986
3. Sacks DB et al., Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. Clin Chem 2002; 48:436-472
4. Hoelzel W et al., IFCC reference system for measurement of haemoglobin A1c in human blood and the national standardization schemes in the United States, Japan and Sweden: a method comparison study. Clin Chem 2004; 50:166-174
5. Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1c Assay into estimated average glucose values. Diabetes Care 2008; 31:1473-1478
6. <http://www.ngsp.org/docs/methods.pdf>
7. <http://www.ngsp.org/protocol.asp>
8. Standards of Medical Care in Diabetes. Recommendations of the American Diabetes Association. Diabetes Care Vol. 34. Jan 2011; 1:S62-S69
9. Parafilm® M is a trademark of Pechiney Plastics Packaging

SYMBOL	MEANING
	Consult Instructions For Use / User Manual
	In vitro Diagnostic Medical Device
	Batch or Lot Number
	Calibration Barcode
	Product Reference
	Control
	Storage Temperature
	Contents Sufficient for "Σ" Number of Tests
	Expiry Date
	Manufacturer
	Caution
	Conformity to the European Directive 98/79/EC on in vitro Diagnostic Medical Devices

## Test Cartridge Storage

### Refrigerated 2 °C to 8 °C (36 °F to 46 °F)

Pouched Quo-Test A1C Test Cartridges are stable until the expiry date only when stored refrigerated. The expiry date is the last day of the month stated on the pouch and kit carton.

### Room Temperature Storage 15 °C to 25 °C (59 °C to 77 °F)

Pouched Quo-Test A1C Test Cartridges can be stored at room temperature for up to 30 days, if still within the expiry date stated on the pouch or carton. Note the date of removal from the refrigerator on the pouch or kit carton. Avoid exposure to direct sunlight.

## Test Cartridge Preparation and Integrity Check

Allow at least 40 minutes for unopened pouches to equilibrate to the analyzer's environmental temperature before use, even if previously stored at room temperature. **Do not** use heaters to warm the cartridges.

Upon opening the pouch, check the color of the silica gel in the sachet against the color chart on the inside back cover of this IFU and explanations below. If the silica gel is completely green, as indicated in figure 4, the cartridge **must not** be used. Contact your local distributor or contact customer support on +49 39203 511 414 or support@ekf-diagnostic.de.

After opening the pouch, the cartridge must be used within 5 minutes

**Please refer to Check Test Cartridge Integrity images 1-4 inside back cover.**

### Figure 1

Silica gel in sachet is orange.  
Cartridge is suitable for use.

### Figure 2

Silica gel in sachet is partially green.  
Cartridge is suitable for use.

### Figure 3

Silica gel in sachet is mostly green.  
Cartridge is suitable for use.

### Figure 4

Silica gel in sachet is completely green.  
**Cartridge is not suitable for use.**

## Intended Use

The Quo-Test A1C Test is intended for the in vitro quantitative determination of glycated hemoglobin (HbA1c) in whole blood obtained from finger prick or from venous blood samples into EDTA tubes.

The Quo-Test A1C Test is indicated in the management and treatment of diabetes mellitus and for monitoring long term glycemic control in patients diagnosed with diabetes. The Quo-Test Analyzer and A1C Test are designed for professional use only.

## Summary and Explanation of the Quo-Test A1C Test

Human red blood cells (erythrocytes) are freely permeable to glucose present in the surrounding liquid (plasma) of the blood.

During the lifetime of the erythrocytes (normally up to 120 days) exposure to blood glucose results in binding of glucose to the hemoglobin A molecule present in erythrocytes. This is referred to as glycated hemoglobin<sup>(1)</sup>.

Chronic elevated blood glucose will result in tissue and organ damage<sup>(2)</sup>.

Good control of blood glucose, evident as lower HbA1c values, has been proven to result in delayed onset and slower progression of these complications<sup>(2,3)</sup>.

## Principle of the Test

The Quo-Test A1C Test Cartridge contains all of the reagents required for the determination of glycated hemoglobin in whole blood.

The Quo-Test A1C Test combines the chemical binding of boronate to glycated hemoglobin with the fluorescent quenching effect that this causes on a fluorescent marker bound to the boronate molecule.

The sample is collected using the Blood Collector which is subsequently placed into the Quo-Test A1C Test Cartridge.

The Test Cartridge is then placed in the Quo-Test Analyzer and the door closed.

The analyzer automatically pushes a mixing bead containing the fluorescent boronate conjugate into the buffer contained within the cartridge and reconstitutes it.

The analyzer then pushes the sample contained on the Blood Collector into the buffer and mixes it to release the hemoglobin from the erythrocytes.

The total hemoglobin concentration is determined from the initial decrease in the fluorescent signal.

The fluorescent boronate conjugate binds to the glycated hemoglobin which is measured by monitoring a decrease in the fluorescence of the active ingredient.

The ratio of glycated hemoglobin to total hemoglobin is determined and the result is presented in up to two user selectable units: % DCCT, mmol/mol IFCC, % JDS, eAG mg/dl or eAG mmol/l.

Linear relationships have been established between these reportable units<sup>(4,5)</sup>:

$$\text{mmol/mol IFCC} = (\% \text{ DCCT} - 2.15) \times 10.929$$

$$\% \text{ JDS} = (0.09274 \times \text{mmol/mol IFCC}) + 1.724$$

$$\text{eAG mg/dl} = (28.7 \times \% \text{ DCCT}) - 46.7$$

$$\text{eAG mmol/l} = (1.59 \times \% \text{ DCCT}) - 2.59$$

eAG values are based on a correlation study linking % DCCT to the patient's average blood glucose concentration, resulting in published formula to derive the eAG<sup>(6)</sup>.

## Cartridge Pack Contents

15 Quo-Test A1C Test Cartridges packed individually into foil pouches an indicating silica gel sachet.

A minimum of 15 Quo-Test Blood Collectors.

1 Calibration Barcode (on inner flap of carton, indicated by the CAL symbol).

1 Instruction For Use leaflet.

## Accessories Required (not included)

Quo-Test Analyzer

Quo-Test A1C Control Kit

Standard Blood Collection Equipment

## Active Ingredients

Fluorescent boronate conjugate

Plastic vial containing ammonium chloride buffer, lysis agent and preservative

Blood Collector treated with EDTA and surfactant

## Calibration of the Quo-Test A1C Test

The Quo-Test A1C System is certified by the NGSP and IFCC<sup>(6)</sup>.

The Quo- Test Analyzer and A1C Test Cartridges have been calibrated using samples provided by the ERL via the NGSP network.

Results obtained using the Quo-Test A1C System are traceable to the IFCC reference method<sup>(7)</sup>.

## Limitations of Use

For hemoglobin concentrations outside of the range of 6.5 g/dL to 20.4 g/dL (4.0 to 12.7 mmol/l) or for HbA1c values outside the range of 4.0 to 15.0 % DCCT (20 to 140 mmol/mol IFCC), a result will not be reported and an error message will be displayed. Consult the Quo-Test Analyzer System User Manual for further details.

eAG values may differ significantly from a patient's glucose concentration if checked at the same time.

As with all HbA1c analyses, the results obtained are dependent upon the patient's clinical history and having red blood cells with a typical lifespan.

## Warnings and Precautions

For *in vitro* diagnostic use only.

**Do not** use Test Cartridges after the expiry date or if the Test Cartridges have not been stored in accordance with the recommendations.

**Do not** use the Test Cartridge if the foil pouch, Test Cartridge or silica gel sachet have been damaged or if the silica gel in the sachet is completely green.

The Test Cartridge contains sodium azide as a preservative. The concentration is < 0.1 %, which is below the level considered to be hazardous in normal use. In case of leakage, avoid contact with eyes and skin. In case of contact, wash with plenty of water.

**Do not** re-use any part of the Test Cartridge.

**Use Gloves** – Used Test Cartridges, sampling equipment, patient samples and controls are potentially infectious. Used Test Cartridges should be disposed of immediately as clinical waste in accordance with local / national regulations.

Operate the Quo-Test Analyzer System according to the User Manual.

## Sample Preparation and Collection

### Finger Prick Blood Sample

The patient's finger must be warm, dry and clean, including being free of substances such as hand cream.

Use a single use lancet (not provided) to prick the finger.

Allow a droplet of blood to form, about the width of the Blood Collector, before sampling.

Fill the Blood Collector, avoiding air bubbles and excess sample on the outside of the tip of the Blood Collector. **Do not** wipe excess blood from the collector.

Place the Blood Collector into the Test Cartridge and the Test Cartridge into the Quo-Test Analyzer within 1 minute.

Operate the Quo-Test Analyzer System according to the User Manual.

## Venous Blood Sample

Venous samples must have been collected into EDTA tubes.

Venous samples can be stored refrigerated at 2 °C to 8 °C (36 °F to 46 °F) for up to 10 days before use and can be used without equilibration to room temperature.

Thoroughly mix the sample by inverting the tube 8 to 10 times before collecting a sample.

Place a droplet of blood, approximately 10 µl to 20 µl, onto a non-metallic and non-absorbent surface such as Parafilm® M<sup>(9)</sup>.

**Do not** attempt to sample directly from the tube.

Fill the Blood Collector, avoiding air bubbles and excess sample on the outside of the tip of the Blood Collector. **Do not** wipe excess blood from the collector.

Place the Blood Collector into the Test Cartridge and the Test Cartridge into the Quo-Test Analyzer within 1 minute.

Operate the Quo-Test Analyzer System according to the User Manual.

## Expected Values

Interpret test results with careful consideration to the patient's age, ethnicity, medical history, clinical examinations and other laboratory results, see limitations of use.

Users should refer to national guidance regarding patient reference ranges, an example is as follows<sup>(8)</sup>:

The ADA recommends an HbA1c range of 4 to 5.7 % DCCT (20 to 39 mmol/mol IFCC) for non-diabetic patients.

Patients with values in the range of 5.7 to 6.0 % DCCT (39 to 42 mmol/mol IFCC) are considered to be at risk of developing diabetes and those in the range of 6.0 to 6.5 % DCCT (42 to 48 mmol/mol IFCC) are considered to be at high risk of developing diabetes.

The ADA recommends a goal of < 7 % DCCT (53 mmol/mol IFCC) for effective management of diabetes and to minimize long term complications.

## Quality Control

The Quo-Test A1C Control Kit must be used with the Quo-Test A1C Test Cartridges to ensure that the system is operating correctly.

Only Quo-Test A1C Controls should be used with the Quo-Test A1C Test.

Please read the IFU included with each Control Kit to ensure that the correct procedure is followed.

Quality control should be run:

- With each new lot of Test Cartridges.
- With each new shipment of Test Cartridges.
- At any time there is concern that the result may be wrong, the Test Cartridges may not have been stored correctly, the user is unfamiliar with the test procedure or there is a concern that a user may be performing the test incorrectly.

Use the Quo-Test A1C Control Kit in accordance with the quality standards set by your local organization or laboratory.

If the Quo-Test A1C Controls do not give the results within the range stated in the control kit IFU, contact your local distributor or contact customer support on +49 39203 511 414 or support@ekf-diagnostic.de

## Performance Characteristics

### Method Comparison

Quo-Test A1C results were compared with two commercially available assays on two separate occasions. Comparison with a POC system gave a linear regression of  $r = 0.971$  ( $y = 0.911 + 0.688x$ ,  $n = 50$ ) and  $r = 0.98$  ( $y = 0.925 + 0.505x$ ,  $n = 48$ ). Comparison with a laboratory HPLC system gave a linear regression of  $r = 0.984$  ( $y = 0.912 + 0.434x$ ,  $n = 50$ ) and  $r = 0.989$  ( $y = 0.977 + 0.08x$ ,  $n = 48$ ).

### Precision

Blood samples with high (9.9 %) and low (5.3 %) A1C values were tested in duplicate every morning and afternoon for twenty days on a single Quo-Test Analyzer. The total imprecision for each level was  $< 3$  % CV respectively.

### Linearity and Test Range

High and low pooled blood samples were used to create eleven linearity standards, which were tested on the Quo-Test A1C Test to establish the linear range. The Quo-Test A1C Test has a linear working range of 4 to 15 % DCCT (20 mmol/mol to 140 mmol/mol IFCC). If a result is outside of this range then the analyzer will display an error message. Refer to the Quo-Test Analyzer System User Manual for further details.

## Interfering Substances

No affect from bilirubin (20 mg/dl), cholesterol (500 mg/dl), creatinine (30 mg/dl), triglycerides (1,500mg/dl), uric acid (20 mg/dl), acetaminophen (20 mg/dl), ascorbic acid (3 mg/dl), dopamine (13 mg/dl), ibuprofen (40 mg/dl), salicylate (50 mg/dl), tetracycline (4mg/dl), tolazamide (100mg/dl) and tolbutamide (100 mg/dl) was found.

## Hemoglobin Variants

The Quo-Test A1C Test was found to be unaffected by labile glycosylated hemoglobin and the following hemoglobin variants: Hb AS, Hb AC, Hb AD, Hb AE, Hb AJ, Hb SS, Hb CC, Hb SC, and elevated fetal hemoglobin (up to 6 %).

It is possible that other substances and/or factors not listed above may interfere and cause false results.

## Glossary of Terms

<b>EDTA</b>	Ethylenediaminetetraacetic Acid
<b>ADA</b>	American Diabetes Association
<b>DCCT</b>	Diabetes Control & Complications Trial
<b>NGSP</b>	National Glycohemoglobin Standardization Program
<b>IFCC</b>	International Federation of Clinical Chemistry
<b>IFU</b>	Instructions For Use
<b>eAG</b>	Estimated Average Glucose
<b>POC</b>	Point Of Care
<b>HPLC</b>	High Pressure Liquid Chromatography

# How to Run a Quo-Test A1C Test

Refer to images 1 - 9 on the inside of the back cover.

## Step 1

Ensure that the analyzer is ready to run a test. Refer to the Quo-Test Analyzer System User Manual for details.

## Step 2

Before using a new lot of Test Cartridges, scan the calibration barcode that is printed on the inner flap of the carton, as indicated by the **CAL** symbol.

Place the barcode scanner over the calibration barcode and press the barcode scanner button until you hear a beep. The analyzer will confirm that the data has been accepted, along with displaying the lot number and expiry date.

## Step 3

Remove a pouched cartridge from the carton and place it next to the analyzer. Allow at least 40 minutes for it to equilibrate to the analyzer's environmental temperature before use, even if previously stored at room temperature.

When ready to perform the analysis, carefully open the pouch and remove the Test Cartridge and silica gel sachet. Check the indicator colour of the silica gel (see colour chart inside back cover).

**Do not** handle the lower part of the cartridge containing the liquid.

**Do not** use the cartridge if condensation is present.

After opening the pouch, the cartridge must be used within 5 minutes.

## Step 4

For finger prick sample collection, the patient's finger must be warm, dry and clean, including being free of substances such as hand cream. Use a single use lancet (not provided) on the patient's finger to get a droplet of blood about the same width as the Blood Collector.

For venous samples previously collected into EDTA tubes, ensure that the sample is thoroughly mixed and use standard sampling equipment to place a droplet of blood, approximately 10 µl to 20 µl, onto a non-metallic and non-absorbent surface such as Parafilm® M<sup>(9)</sup>. **Do not** attempt to sample directly from the tube.

## Step 5

Lightly touch the pointed tip of the Blood Collector to the top surface of the droplet of blood, as per the image. The blood will be taken up by the collector. The blood must completely fill the slot.

**Note:** Avoid taking up air bubbles and excess sample on the outside of the tip of the Blood Collector. **Do not** wipe excess blood from the collector.

## Step 6

Cartridge, making sure that the top of the Blood Collector is flush with the top of the Test Cartridge and not sticking out. **Do not** push the Blood Collector into the Test Cartridge, as this may cause an error message to be displayed instead of a result. Refer to the Quo-Test Analyzer System User Manual for further details.

## Step 7

After the Test Cartridge has been prepared, open the door of the analyzer and place the Test Cartridge into the test chamber, Making sure that the Test Cartridge is firmly seated.

The Test Cartridge must be inserted into the analyzer within 1 minute of the blood sample being collected.

## Step 8

Pull the red slide towards the front of the analyzer and then close the door. The test will start automatically. You will be prompted to scan the Operator and Patient ID if these features have been enabled.

**Note:** If the analyzer displays 'Scan Lot \*\*\*\*\* Data Card' then scan the calibration barcode (see Step 2).

## Step 9

When the test has finished the result will be displayed on the screen. Open the door and push the red slide towards the rear of the analyzer. Remove the Test Cartridge and dispose of it as clinical waste. Upon closing the analyzer door, the display will return to the home screen and the analyzer will be ready to run another sample.

**Note:** Refer to the Quo-Test Analyzer User Manual to access previous results in memory.

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## **Ordliste**

<b>EDTA</b>	Etylendiamintetraeddiksyre
<b>ADA</b>	American Diabetes Association
<b>DCCT</b>	Diabetes Control and Complications Trial
<b>NGSP</b>	National Glycohemoglobin Standardization Program
<b>IFCC</b>	International Federation of Clinical Chemistry
<b>IFU</b>	Instructions For Use (bruksanvisning)
<b>eAG</b>	Estimated Average Glucose (anslått gjennomsnittlig glukose)
<b>POC</b>	Point Of Care
<b>HPLC</b>	High Pressure Liquid Chromatography (væskechromatografi)

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