



Special features and codes tipstim® pulse generator

not explicitly listed in the instruction manual
(as of April 12, 2016)

1. Entering the menu

- Switch on the device by pressing the "ON/OFF"- button.
- After 1-2 seconds the tipstim® logo will appear on the screen.
- Immediately press the "OK"- button and keep it pressed.
- As soon as the menu heading "code" appears on the screen release the "OK"- button.
- The code to enter the menu is always "0392" (reminder: the first 3 digits of the postal code of the manufacturer on the back side of the device headed by "0").

Every code is always entered via the selection buttons "+"- and "-" (either of the left or right channel). The correct entry has to be confirmed by pressing the "OK"- button.

2. Programme selection (P1 or P2)

- Make sure the device is turned off.
- Press "OK" button and "+" - button of the right channel at a time and keep it pressed.
- Now press the "ON/OFF" - button (while keeping "OK" button and "+" button of the right channel pressed).
- As soon as the menu heading "Choose Program..." appears on the screen release all buttons.
- Select the desired program by pressing "+" and "-" – buttons of the left or right channel (Program 1: standard program tipstim®, Program 2: special program to recede spasticity).
- Confirm selected program by pressing "OK" button.

3. **One-time deletion of the “glove counter”** (following messages after 90 – 100 applications of tipstim® glove - to be used by the patient when changing tipstim® glove for a new one)

**Please change
the Glove
Only 2 Sessions
remaining** **ATTENTION!
Please change
glove.**

After changing tipstim® glove for a new one messages can be deleted as follows:

- Make sure the device is turned off.
- Press “-“ - button of the left channel and keep it pressed.
- Now press the “ON/OFF“ - button (while keeping “-“ - button of the left channel pressed).
- As soon as the menu heading “code” appears on the screen release all buttons.
- Enter code “0314”.
- Confirm question “RESET (Y/N)” by pressing the “OK” - button.
- Counter is reset to 0.

This code won't be provided to the patient unless he gets a fresh tipstim® glove after 100 applications.

4. **Permanent deletion of the “glove counter” mentioned under 3.** (to be used by doctors/therapists)

- Make sure the device is turned off.
- Press “-“ - buttons of both channels and “OK” - button at a time and keep it pressed (i.e. 3 buttons at a time)
- Now press the “ON/OFF“ - button (while keeping “-“ - buttons of both channels and “OK” – button pressed).
- As soon as the menu heading “code” appears on the screen release all buttons.
- Enter code “0315”.
- “Glove Counter” can be turned on/off by using “+“ and “-“ - buttons of the right channel.

5. Complete reset to factory settings

- Make sure the device is turned off.
- Press “+” and „-“ - buttons of the right channel at a time and keep it pressed.
- Now press the “ON/OFF” - button (while keeping “+“- und “-“ - buttons of the right channel pressed).
- As soon as the menu heading “code” appears on the screen release all buttons.
- Enter code “6941”.
- Confirm question “RESET (Y/N)” by pressing the “OK” - button.
- The correct reset is confirmed by 2 confirmation beeps.

6. Checking the remaining battery power

- Make sure the device is turned off.
- Press “+” - button of the right channel and keep it pressed.
- Now press the “ON/OFF” - button once(while keeping “+“ - button of the right channel pressed at a time).
- As long as you keep the “+” - button of the right channel pressed display will show you the remaining battery power as a bar graph.