

Quo-Test™ HbA1C Analyser	Serial number: SN01
---------------------------------	----------------------------

System Quality Control Log sheet - NB: For each Quo-Test™ HbA1C analyser you must use a separate log sheet

Both controls, **Level 1 Normal** (L1; white (w)) and **Level 2 Abnormal** (L2; black (b)) should be used for a Quality Control test

Date	Reason for test (routine; new shipment or lot of test cartridges; suspected error; other - please define)	Lot number of control solutions	Expiry date of control solutions	Date of opening	Batch number of test cartridges	Expiry date of test cartridges	Expected results range				Test result				Pass (P)/ Fail (F)	Operator initials
							L1 (w)		L2 (b)		L1 (w)		L2 (b)			
							DCCT (%)	IFCC (mmol/mol)	DCCT (%)	IFCC (mmol/mol)	DCCT (%)	IFCC (mmol/mol)	DCCT (%)	IFCC (mmol/mol)		
09/11/12	New shipment of test cartridges	CS12b	2016-11	9/11/15	020089	2016-12	6.0 – 7.7	42.1 – 60.7	10.3 – 12.6	89.1 – 114.2	6.9	51.9	11.5	102.2	P	JAC